

Product Datasheet

Human Regenerating Islet Derived Protein 3 Gamma (REG3g) ELISA Kit (orb1204899)

Description	Human Regenerating Islet Derived Protein 3 Gamma (REG3g) ELISA Kit
Reactivity	Human
Range	31.2-2000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q6UW15
Sensitivity	14.5pg/mL
Entrez	130120
Expiration Date	Please enquire.