

Product Datasheet

Human Ras Association Domain Containing Protein 1 (RASSF1) ELISA Kit (orb1204905)

Description	Human Ras Association Domain Containing Protein 1 (RASSF1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q9NS23
Sensitivity	0.055ng/mL
Entrez	11186
Expiration Date	Please enquire.