

Product Datasheet

Human Retinoic Acid Receptor Responder 1 (RARRES1) ELISA Kit (orb1204907)

Description	Human Retinoic Acid Receptor Responder 1 (RARRES1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P49788
Sensitivity	0.066ng/mL
Entrez	5918
Expiration Date	Please enquire.