

Product Datasheet

Human Protein Tyrosine Phosphatase Like Protein A (PTPLA) ELISA Kit (orb1204932)

Description	Human Protein Tyrosine Phosphatase Like Protein A (PTPLA) ELISA Kit
Reactivity	Human
Range	0.781-50ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	BOYJ81
Sensitivity	0.25ng/mL
Entrez	9200
Expiration Date	Please enquire.