

Product Datasheet

Human Protein Tyrosine Phosphatase Type IVA 3 (PTP4A3) ELISA Kit (orb1204933)

Description	Human Protein Tyrosine Phosphatase Type IVA 3 (PTP4A3) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	075365
Sensitivity	0.118ng/mL
Entrez	11156
Expiration Date	Please enquire.