

Product Datasheet

Human Proteasome Subunit Beta Type 6 (PSMb6) ELISA Kit (orb1204941)

Description	Human Proteasome Subunit Beta Type 6 (PSMb6) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P28072
Sensitivity	0.053ng/mL
Entrez	5694
Expiration Date	Please enquire.