

Product Datasheet

Human Prosaposin (PSAP) ELISA Kit (orb1204949)

Description	Human Prosaposin (PSAP) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cerebrospinal fluid, seminal plasma, cell culture supernates or other biological fluids.
Uniprot ID	P07602
Sensitivity	0.055ng/mL
Entrez	5660
Expiration Date	Please enquire.