

Product Datasheet

Human Peroxiredoxin 5 (PRDX5) ELISA Kit (orb1204961)

Description	Human Peroxiredoxin 5 (PRDX5) ELISA Kit
Reactivity	Human
Range	0.781-50ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cerebrospinal fluid or other biological fluids.
Uniprot ID	P30044
Sensitivity	0.34ng/mL
Entrez	25824
Expiration Date	Please enquire.