

## Product Datasheet

# Human Serine/Threonine-Protein Phosphatase 2A 56 kDa Regulatory Subunit Alpha Isoform (PP2Aa) ELISA Kit (orb1204977)

|                          |                                                                                                            |
|--------------------------|------------------------------------------------------------------------------------------------------------|
| <b>Description</b>       | Human Serine/Threonine-Protein Phosphatase 2A 56 kDa Regulatory Subunit Alpha Isoform (PP2Aa) ELISA<br>... |
| <b>Reactivity</b>        | Human                                                                                                      |
| <b>Range</b>             | 78.1-5000pg/mL                                                                                             |
| <b>Note</b>              | For research use only                                                                                      |
| <b>Application notes</b> | sandwich                                                                                                   |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates<br>or other biological fluids.                                           |
| <b>Uniprot ID</b>        | <b>P30153</b>                                                                                              |
| <b>Sensitivity</b>       | 30.9pg/mL                                                                                                  |
| <b>Entrez</b>            | <b>5518</b>                                                                                                |
| <b>Expiration Date</b>   | Please enquire.                                                                                            |