

Product Datasheet

Human Pim-1 Oncogene (PIM1) ELISA Kit (orb1205016)

Description	Human Pim-1 Oncogene (PIM1) ELISA Kit
Reactivity	Human
Range	0.781-50ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P11309
Sensitivity	0.25ng/mL
Entrez	5292
Expiration Date	Please enquire.