

Product Datasheet

Human Peptidoglycan Recognition Protein 1 (PGLYRP1) ELISA Kit (orb1205023)

Description	Human Peptidoglycan Recognition Protein 1 (PGLYRP1) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, cell lysates or other biological fluids.
Uniprot ID	075594
Sensitivity	0.115ng/mL
Entrez	8993
Expiration Date	Please enquire.