

Product Datasheet

Human Phosphatidylethanolamine Binding Protein 1 (PEBP1) ELISA Kit (orb1205030)

Description	Human Phosphatidylethanolamine Binding Protein 1 (PEBP1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P30086
Sensitivity	0.058ng/mL
Entrez	5037
Expiration Date	Please enquire.