

Product Datasheet

Human Prolyl-4-Hydroxylase Alpha Polypeptide III (P4Ha3) ELISA Kit (orb1205071)

Description	Human Prolyl-4-Hydroxylase Alpha Polypeptide III (P4Ha3) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q7Z4N8
Sensitivity	0.066ng/mL
Entrez	283208
Expiration Date	Please enquire.