

## Product Datasheet

# Human Prolyl-4-Hydroxylase Alpha Polypeptide I (P4Ha1) ELISA Kit (orb1205073)

|                          |                                                                        |
|--------------------------|------------------------------------------------------------------------|
| <b>Description</b>       | Human Prolyl-4-Hydroxylase<br>Alpha Polypeptide I (P4Ha1)<br>ELISA Kit |
| <b>Reactivity</b>        | Human                                                                  |
| <b>Range</b>             | 0.312-20ng/mL                                                          |
| <b>Note</b>              | For research use only                                                  |
| <b>Application notes</b> | sandwich                                                               |
| <b>Sample Types</b>      | tissue homogenates or other<br>biological fluids.                      |
| <b>Uniprot ID</b>        | <b>P13674</b>                                                          |
| <b>Sensitivity</b>       | 0.112ng/mL                                                             |
| <b>Entrez</b>            | <b>5033</b>                                                            |
| <b>Expiration Date</b>   | Please enquire.                                                        |