

## Product Datasheet

# Human Oxidative Stress Induced Growth Inhibitor 1 (OSGIN1) ELISA Kit (orb1205082)

|                          |                                                                      |
|--------------------------|----------------------------------------------------------------------|
| <b>Description</b>       | Human Oxidative Stress Induced Growth Inhibitor 1 (OSGIN1) ELISA Kit |
| <b>Reactivity</b>        | Human                                                                |
| <b>Range</b>             | 0.156-10ng/mL                                                        |
| <b>Note</b>              | For research use only                                                |
| <b>Application notes</b> | sandwich                                                             |
| <b>Sample Types</b>      | tissue homogenates, cell lysates or other biological fluids.         |
| <b>Uniprot ID</b>        | <a href="#">Q9UJX0</a>                                               |
| <b>Sensitivity</b>       | 0.064ng/mL                                                           |
| <b>Entrez</b>            | <a href="#">29948</a>                                                |
| <b>Expiration Date</b>   | Please enquire.                                                      |