

Product Datasheet

Human Osteoclast Associated Receptor (OSCAR) ELISA Kit (orb1205084)

Description	Human Osteoclast Associated Receptor (OSCAR) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q8IYS5
Sensitivity	0.136ng/mL
Entrez	126014
Expiration Date	Please enquire.