

Product Datasheet

Human Oxysterol Binding Protein Like Protein 8 (OSBPL8) ELISA Kit (orb1205085)

Description	Human Oxysterol Binding Protein Like Protein 8 (OSBPL8) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q9BZF1
Sensitivity	0.138ng/mL
Entrez	114882
Expiration Date	Please enquire.