

Product Datasheet

Human Oxoglutarate Dehydrogenase (OGDH) ELISA Kit (orb1205090)

Description	Human Oxoglutarate Dehydrogenase (OGDH) ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q02218
Sensitivity	0.59ng/mL
Entrez	4967
Expiration Date	Please enquire.