

Product Datasheet

Human 2',5'-Oligoadenylate Synthetase Like Protein (OASL) ELISA Kit (orb1205093)

Description	Human 2',5'-Oligoadenylate Synthetase Like Protein (OASL) ELISA Kit
Reactivity	Human
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, cell lysates or other biological fluids.
Uniprot ID	Q15646
Sensitivity	6.2pg/mL
Entrez	8638
Expiration Date	Please enquire.