

## Product Datasheet

# Human 2',5'-Oligoadenylate Synthetase 2 (OAS2) ELISA Kit (orb1205095)

|                          |                                                               |
|--------------------------|---------------------------------------------------------------|
| <b>Description</b>       | Human 2',5'-Oligoadenylate Synthetase 2 (OAS2) ELISA Kit      |
| <b>Reactivity</b>        | Human                                                         |
| <b>Range</b>             | 3.12-200ng/mL                                                 |
| <b>Note</b>              | For research use only                                         |
| <b>Application notes</b> | sandwich                                                      |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates or other biological fluids. |
| <b>Uniprot ID</b>        | <b>P29728</b>                                                 |
| <b>Sensitivity</b>       | 1.24ng/mL                                                     |
| <b>Entrez</b>            | <b>4939</b>                                                   |
| <b>Expiration Date</b>   | Please enquire.                                               |