

Product Datasheet

Human 2',5'-Oligoadenylate Synthetase 1 (OAS1) ELISA Kit (orb1205096)

Description	Human 2',5'-Oligoadenylate Synthetase 1 (OAS1) ELISA Kit
Reactivity	Human
Range	4.68-300ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, cell lysates or other biological fluids.
Uniprot ID	P00973
Sensitivity	1.8ng/mL
Entrez	4938
Expiration Date	Please enquire.