

Product Datasheet

Human Nuclear Receptor Related Protein 1 (NURR1) ELISA Kit (orb1205099)

Description	Human Nuclear Receptor Related Protein 1 (NURR1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P43354
Sensitivity	0.053ng/mL
Entrez	4929
Expiration Date	Please enquire.