

## Product Datasheet

### Human Nucleoporin 62kDa (NUP62) ELISA Kit (orb1205100)

|                          |                                                   |
|--------------------------|---------------------------------------------------|
| <b>Description</b>       | Human Nucleoporin 62kDa<br>(NUP62) ELISA Kit      |
| <b>Reactivity</b>        | Human                                             |
| <b>Range</b>             | 0.312-20ng/mL                                     |
| <b>Note</b>              | For research use only                             |
| <b>Application notes</b> | sandwich                                          |
| <b>Sample Types</b>      | tissue homogenates or<br>other biological fluids. |
| <b>Uniprot ID</b>        | <b>P37198</b>                                     |
| <b>Sensitivity</b>       | 0.112ng/mL                                        |
| <b>Entrez</b>            | <b>23636</b>                                      |
| <b>Expiration Date</b>   | Please enquire.                                   |