

Product Datasheet

Human Neurotensin Receptor 1, High Affinity (NTSR1) ELISA Kit (orb1205105)

Description	Human Neurotensin Receptor 1, High Affinity (NTSR1) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P30989
Sensitivity	0.126ng/mL
Entrez	4923
Expiration Date	Please enquire.