

Product Datasheet

Human Neuregulin 3 (NRG3) ELISA Kit (orb1205115)

Description	Human Neuregulin 3 (NRG3) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P56975
Sensitivity	0.063ng/mL
Entrez	10718
Expiration Date	Please enquire.