

Product Datasheet

Human Nuclear Respiratory Factor 1 (NRF1) ELISA Kit (orb1205117)

Description	Human Nuclear Respiratory Factor 1 (NRF1) ELISA Kit
Reactivity	Human
Range	78.1-5000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q16656
Sensitivity	30.98pg/mL
Entrez	4899
Expiration Date	Please enquire.