

Product Datasheet

Human Neuropeptide Y Receptor Y5 (NPY5R) ELISA Kit (orb1205120)

Description	Human Neuropeptide Y Receptor Y5 (NPY5R) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q15761
Sensitivity	0.108ng/mL
Entrez	4889
Expiration Date	Please enquire.