

Product Datasheet

Human Neuropeptide Y Receptor Y2 (NPY2R) ELISA Kit (orb1205121)

Description	Human Neuropeptide Y Receptor Y2 (NPY2R) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P49146
Sensitivity	0.114ng/mL
Entrez	4887
Expiration Date	Please enquire.