

## Product Datasheet

### Human Nucleophosmin (NPM) ELISA Kit (orb1205130)

|                          |                                                                 |
|--------------------------|-----------------------------------------------------------------|
| <b>Description</b>       | Human Nucleophosmin (NPM)<br>ELISA Kit                          |
| <b>Reactivity</b>        | Human                                                           |
| <b>Range</b>             | 0.312-20ng/mL                                                   |
| <b>Note</b>              | For research use only                                           |
| <b>Application notes</b> | sandwich                                                        |
| <b>Sample Types</b>      | tissue homogenates, cell lysates<br>or other biological fluids. |
| <b>Uniprot ID</b>        | <b>P06748</b>                                                   |
| <b>Sensitivity</b>       | 0.128ng/mL                                                      |
| <b>Entrez</b>            | <b>4869</b>                                                     |
| <b>Expiration Date</b>   | Please enquire.                                                 |