

Product Datasheet

Human Non Metastatic Cells 1, Protein NM23A Expressed In (NME1) ELISA Kit (orb1205136)

Description	Human Non Metastatic Cells 1, Protein NM23A Expressed In (NME1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	P15531
Sensitivity	0.059ng/mL
Entrez	4830
Expiration Date	Please enquire.