

Product Datasheet

Human Nucleosome Assembly Protein 1 Like Protein 1 (NAP1L1) ELISA Kit (orb1205172)

Description	Human Nucleosome Assembly Protein 1 Like Protein 1 (NAP1L1) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P55209
Sensitivity	0.093ng/mL
Entrez	4673
Expiration Date	Please enquire.