

Product Datasheet

Human Myosin Heavy Chain 7, Cardiac Muscle, Beta (MYH7) ELISA Kit (orb1205187)

Description	Human Myosin Heavy Chain 7, Cardiac Muscle, Beta (MYH7) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P12883
Sensitivity	0.053ng/mL
Entrez	4625
Expiration Date	Please enquire.