

Product Datasheet

Human Myosin Heavy Chain 4, Skeletal Muscle (MYH4) ELISA Kit (orb1205189)

Description	Human Myosin Heavy Chain 4, Skeletal Muscle (MYH4) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q9Y623
Sensitivity	0.127ng/mL
Entrez	4622
Expiration Date	Please enquire.