

Product Datasheet

Human Maternal Embryonic Leucine Zipper Kinase (MELK) ELISA Kit (orb1205236)

Description	Human Maternal Embryonic Leucine Zipper Kinase (MELK) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q14680
Sensitivity	0.105ng/mL
Entrez	9833
Expiration Date	Please enquire.