

Product Datasheet

Human Leucine Rich Repeat In FLII Interacting Protein 1 (LRRFIP1) ELISA Kit (orb1205267)

Description	Human Leucine Rich Repeat In FLII Interacting Protein 1 (LRRFIP1) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q32MZ4
Sensitivity	0.127ng/mL
Entrez	9208
Expiration Date	Please enquire.