

Product Datasheet

Human Lysophosphatidic Acid Receptor 3 (LPAR3) ELISA Kit (orb1205277)

Description	Human Lysophosphatidic Acid Receptor 3 (LPAR3) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q9UBY5
Sensitivity	0.105ng/mL
Entrez	23566
Expiration Date	Please enquire.