

Product Datasheet

Human Lysyl Oxidase Like Protein 4 (LOXL4) ELISA Kit (orb1205279)

Description	Human Lysyl Oxidase Like Protein 4 (LOXL4) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q96JB6
Sensitivity	0.053ng/mL
Entrez	84171
Expiration Date	Please enquire.