

Product Datasheet

Human Lemur Tyrosine Kinase 3 (LMTK3) ELISA Kit (orb1205286)

Description	Human Lemur Tyrosine Kinase 3 (LMTK3) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q96Q04
Sensitivity	0.054ng/mL
Entrez	114783
Expiration Date	Please enquire.