

Product Datasheet

Human Leucine Rich Repeat LGI Family, Member 3 (LGI3) ELISA Kit (orb1205295)

Description	Human Leucine Rich Repeat LGI Family, Member 3 (LGI3) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	Q8N145
Sensitivity	0.055ng/mL
Entrez	203190
Expiration Date	Please enquire.