

Product Datasheet

Human Leprecan Like Protein 1 (LEPREL1) ELISA Kit (orb1205296)

Description	Human Leprecan Like Protein 1 (LEPREL1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q8IVL5
Sensitivity	0.063ng/mL
Entrez	55214
Expiration Date	Please enquire.