

Product Datasheet

Human Inter Alpha-Globulin Inhibitor H2 (ITIH2) ELISA Kit (orb1205355)

Description	Human Inter Alpha-Globulin Inhibitor H2 (ITIH2) ELISA Kit
Reactivity	Human
Range	78.1-5000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P19823
Sensitivity	35pg/mL
Entrez	3698
Expiration Date	Please enquire.