

Product Datasheet

Human Inosine 5'-Monophosphate Dehydrogenase 2 (IMPDH2) ELISA Kit (orb1205368)

Description	Human Inosine 5'-Monophosphate Dehydrogenase 2 (IMPDH2) ELISA Kit
Reactivity	Human
Range	0.781-50ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P12268
Sensitivity	0.27ng/mL
Entrez	3615
Expiration Date	Please enquire.