

Product Datasheet

Human 5-Hydroxytryptamine Receptor 2A (HTR2A) ELISA Kit (orb1205382)

Description	Human 5-Hydroxytryptamine Receptor 2A (HTR2A) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P28223
Sensitivity	0.057ng/mL
Entrez	3356
Expiration Date	Please enquire.