

Product Datasheet

Human Heat Shock Protein 90kDa Alpha B1 (HSP90aB1) ELISA Kit (orb1205394)

Description	Human Heat Shock Protein 90kDa Alpha B1 (HSP90aB1) ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P08238
Sensitivity	0.67ng/mL
Entrez	3326
Expiration Date	Please enquire.