

Product Datasheet

Human Histidine Rich Calcium Binding Protein (HRC) ELISA Kit (orb1205409)

Description	Human Histidine Rich Calcium Binding Protein (HRC) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P23327
Sensitivity	0.123ng/mL
Entrez	3270
Expiration Date	Please enquire.