

Product Datasheet

Human High Mobility Group AT Hook Protein 1 (HMGA1) ELISA Kit (orb1205417)

Description	Human High Mobility Group AT Hook Protein 1 (HMGA1) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P17096
Sensitivity	0.124ng/mL
Entrez	3159
Expiration Date	Please enquire.