

Product Datasheet

Human Guanylate Cyclase Activator 2A (GUCA2A) ELISA Kit (orb1205445)

Description	Human Guanylate Cyclase Activator 2A (GUCA2A) ELISA Kit
Reactivity	Human
Range	31.2-2000fmol/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q02747
Sensitivity	8.78fmol/mL
Entrez	2980
Expiration Date	Please enquire.