

## Product Datasheet

# Human Glycosylphosphatidylinositol Anchored Molecule Like Protein (GML) ELISA Kit (orb1205475)

|                          |                                                                                   |
|--------------------------|-----------------------------------------------------------------------------------|
| <b>Description</b>       | Human Glycosylphosphatidylinositol Anchored Molecule Like Protein (GML) ELISA Kit |
| <b>Reactivity</b>        | Human                                                                             |
| <b>Range</b>             | 0.156-10ng/mL                                                                     |
| <b>Note</b>              | For research use only                                                             |
| <b>Application notes</b> | sandwich                                                                          |
| <b>Sample Types</b>      | tissue homogenates or other biological fluids.                                    |
| <b>Uniprot ID</b>        | <a href="#">Q99445</a>                                                            |
| <b>Sensitivity</b>       | 0.055ng/mL                                                                        |
| <b>Entrez</b>            | <a href="#">2767</a>                                                              |
| <b>Expiration Date</b>   | Please enquire.                                                                   |