

Product Datasheet

Human Gamma-Glutamyl Hydrolase (gGH) ELISA Kit (orb1205495)

Description	Human Gamma-Glutamyl Hydrolase (gGH) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q92820
Sensitivity	0.056ng/mL
Entrez	8836
Expiration Date	Please enquire.