

Product Datasheet

Human Growth Differentiation Factor 3 (GDF3) ELISA Kit (orb1205497)

Description	Human Growth Differentiation Factor 3 (GDF3) ELISA Kit
Reactivity	Human
Range	78.1-5000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q9NR23
Sensitivity	34pg/mL
Entrez	9573
Expiration Date	Please enquire.